



HOTPLATE BER1BY1

HOTPLATE BER1BY1

NEW HOT PLATE, GRAPHITE HOT PLATE



1. hot plate or Graphite hot plate, uniform surface temperature.
2. up fast and uniform.
3. Digital hot plate: digital display indicating temperature of hot plate clearly.

SPECIFICATIONS

Model	BER1BY1
Function	
Max Temp. (°C)	320°C
Structure	
Outer Shell	Iron plate welding, surface spraying
Hot plate	Aluminum tensile
Heating Power (W)	800
Setting Method	Scale
Timer	No / Digital
Temperature Sensor	K model
Temp. Accuracy	±10%
Plate Temp. Uniformity	±15%
Specification	
Exterior Size (W*L*H) (mm)	340*370*220
Packing Size (W*L*H) (mm)	520*380*250
Plate Size (mm)	300*300
Plate Load	20KG
Power Supply (50/60Hz)	220V
NW/GW (kg)	9/10
Alt Name	New Hot Plate





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com