



WATER BATH BER1CB1

WATER BATH BER1CB1

TRANSPARENT WATER BATH



1. Using innovative design structure, effectively prevent steam leakage.
2. Intelligent temperature control instrument with high brightness digital display, light-touch type key has the advantages of simple operation and high precision temperature control.
3. Shell, interior and cover of the product are all made of transparent engineering plastics, can clearly see the object to be heated.
4. Off when water shortage to ensure the safe operation of the product.

SPECIFICATIONS

Model	BER1CB1
Capacity	2.6L
Voltage (V)	200-240V/100-120V
Power (kW)	0.4
Temp. Motion (°C)	± 0.5
Temp. Range (°C)	Rt+5~100
Temp. Sensitivity (°C)	±± 1
Shelves Height (mm)Distance from bottom	100mm
Shelves Holes Diameter	Φ18.7mm*16
Inner Chamber Size (L*W*H)	143*130*150
Exterior Size (L*W*H)	254*190*193
Packing Size (L*W*H)	295*270*275
N.W.(KG)	2.1
G.W. (KG)	2.7
Alt Name	Transparent Water Bath



APPLICATIONS

Used for precision constant temperature and auxiliary.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com