



WATER BATH BER1CB4

WATER BATH BER1CB4

TRANSPARENT WATER BATH



1. Using innovative design structure, effectively prevent steam leakage.
2. Intelligent temperature control instrument with high brightness digital display, light-touch type key has the advantages of simple operation and high precision temperature control.
3. Shell, interior and cover of the product are all made of transparent engineering plastics, can clearly see the object to be heated.
4. Off when water shortage to ensure the safe operation of the product.

SPECIFICATIONS

| Model | BER1CB4 |
|---|------------------------|
| Capacity | 12L |
| Voltage (V) | 200-240V/100-120V |
| Power (kW) | 1 |
| Temp. Motion (°C) | ± 0.5 |
| Temp. Range (°C) | Rt+5~100 |
| Temp. Sensitivity (°C) | ±± 1 |
| Shelves Height (mm)Distance from bottom | 150mm |
| Shelves Holes Diameter | Φ18.7mm*48 |
| Inner Chamber Size (L*W*H) | 285*225*200 |
| Exterior Size (L*W*H) | 401*276*242 |
| Packing Size (L*W*H) | 450*350*300 |
| N.W.(KG) | 3.9 |
| G.W. (KG) | 4.8 |
| Alt Name | Transparent Water Bath |



APPLICATIONS

Used for precision constant temperature and auxiliary.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com