



PORTABLE STERILIZER BER1CJ2

PORTABLE STERILIZER BER1CJ2

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Double scale indication pressure gauge.
3. Electric heated.
4. Easy to operate, safe and reliable.

SPECIFICATIONS

Model	BER1CJ2
Capacity	24LΦ280*390mm
Working Pressure	0.14-0.16Mpa
Working Temperature	126°C
Max.Safety Pressure	0.165Mpa
Heat Average	$\pm 0.1^{\circ}\text{C}$
Power	2KW/AC200V.50Hz
Packing Size	410*410*630mm
G.W/N.W	19kg/17kg
Alt Name	Portable Pressure Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com