



PORTABLE STERILIZER BER1CK2

PORTABLE STERILIZER BER1CK2

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Electric or LPG heated.
3. Double scale indication pressure gauge.
4. Easy to operate, safe and reliable.

SPECIFICATIONS

Model	BER1CK2
Capacity	24L
Working Pressure	0.14-0.16Mpa
Working Temperature	126°C
Max.Safety Pressure	0.165Mpa
Power	AC220V.50Hz/2Kw
Packing Size(mm)	410*410*550
G.W/N.W	17/15KG
Alt Name	Portable Pressure Steam Sterilizer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com