



## ACCESSORIES BER1CO1

## ACCESSORIES BER1C01

### BASIC POWER SUPPLY



1. Molding shell,touch keys,microprocessor intelligent control.
2. With automatic memory function.
3. With standard,timing operation function.
4. Automatic detection of no-load,over-load.

## SPECIFICATIONS

|                |                              |
|----------------|------------------------------|
| Model          | BER1C01                      |
| Type of Output | Constant-Voltage, Current    |
| Output range   | 5~300V, 1~300mA, Maximum 90W |
| Increment      | 1V, 1mA                      |
| Timer Range    | 1min~9hr, 59min              |
| Display        | LED                          |
| Output jacks   | Two sets of output jacks     |
| Alt Name       | Basic Power Supply           |



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)