



MINI CENTRIFUGE BER1CW1

MINI CENTRIFUGE BER1CW1

DIGITAL CENTRIFUGE



1. Fixed angle rotor with capacity for 6 tubes of 15 mL (dimensions 16x100 mm)
2. Transparent polycarbonate lid for visualization and control of rotor.
3. LED display for adjustment and visualization of operation time and speed.
4. switch that avoids centrifugation if the lid is open and that automatically stops operation in case of opening the lid.
5. suspension motor to assure a soft operation.
6. Sucker feet for a good fixation to the working surface.

SPECIFICATIONS

Model	BER1CW1
Control	Digital
Rotor	Fixed-angle
Speed	1000-4000 rpm (steps of 500 rpm)
Time	1-59 min (steps of 1 min)
Capacity	6x15 ml
Max. RCF	1790 g
Consumption	80 W
Power	220 V - 50 Hz
Alt Name	Digital Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com