



LABORATORY LOW SPEED CENTRIFUGE BER1CX1

LABORATORY LOW SPEED CENTRIFUGE BER1CX1

DIGITAL TABLETOP LOW SPEED CENTRIFUGE



- 1.DC brushless motor drive, digital display speed and running time
- 2.Mechanical door lock
- 3.Used for fat graft purification, PRP and other cosmetic surgery fields.

SPECIFICATIONS

Model	BER1CX1
Maximum speed	4000r/min
Maximum RCF	2350 x g
Max volume	12 x 15ml
Power supply	AC 220V 50HZ
Timer	1 ~ 99min
Noise	≤55dBA
Exterior Size (mm)	270 x 290 x 230
Packing Size (mm)	350 x 340 x 360
N.W / G.W (KG)	8 / 9
Alt Name	Digital Tabletop Low Speed Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com