



WATER BATH BER1DL2

WATER BATH BER1DL2

THERMOSTATIC WATER BATH



Cold-rolling steel spraying exterior.
Stainless steel Inner tank and cover.
Intelligent Temperature control instrument adopts PID control, dual digital display, with functions of timing of operation (0-9999 min), fixed value operation, automatically stop, realization of over-temperature alarm, power-off compensation, power failure memory, error correction and Menu lock.

SPECIFICATIONS

Model	BER1DL2
Voltage	200-240V/100-120V
Power(W)	500
Temmotion (°C)	±0.5
Tem-range (°C)	Room Temp.±5~100°C
Tem. Sensitivity (°C)	≤±1
Display error (°C)	≤±2.5
Chamber Size (mm)	325*168*120
Products Size (mm)	510*240*190
Packing Size (mm)	555*320*230
N.W.(KG)	5
G.W.(KG)	6.5
Alt Name	Thermostatic water bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com