



WATER BATH BER1DL7 (BBWA-107)

WATER BATH BER1DL7

THERMOSTATIC WATER BATH



Cold-rolling steel spraying exterior.
Stainless steel Inner tank and cover.
Intelligent Temperature control instrument adopts PID control, dual digital display, with functions of timing of operation (0-9999 min), fixed value operation, automatically stop, realization of over-temperature alarm, power-off compensation, power failure memory, error correction and Menu lock.

SPECIFICATIONS

Model	BER1DL7
Old Model	BBWA-107
Voltage	200-240V/100-120V
Power(W)	2000
Temmotion (°C)	±0.5
Tem-range (°C)	Room Temp.±5~100°C
Tem. Sensitivity (°C)	≤±1
Display error (°C)	≤±2.5
Chamber Size (mm)	635*325*120
Products Size (mm)	830*380*190
Packing Size (mm)	870*480*235
N.W.(KG)	9
G.W.(KG)	10
Alt Name	Thermostatic water bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com