



LABORATORY HIGH SPEED CENTRIFUGE BET-2703

LABORATORY HIGH SPEED CENTRIFUGE BET-2703

HANDHELD CENTRIFUGES



- Delicate in appearance and small in volume, easy to operate and highly stable.
- Patented mute technology allows the machine to operate at ultra-low noise.
- Micro DC motors and carbon fiber rotors.

SPECIFICATIONS

Model	BET-2703
Max speed	10,000r/min
Max RCF	5,000xg
Rotor configuration	6x1.5/0.5/0.2mL2x8x0.2mL
Dimensions (LxWxH)	170mmx150mmx130mm
Weight	1.25kg
Alt Name	Handheld Centrifuges



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com