



MINI CENTRIFUGE BET-2801

MINI CENTRIFUGE BET-2801

MICROPLATE CENTRIFUGE



- This compact, economical and mini Microplate Centrifuge is an instantaneous centrifuge specially designed for 96-well microplates, which can quickly spin down droplets.
- Most microplate centrifuges on the market are bulky, taking up a lot of lab space.
- This microplate centrifuge has a unique and compact design. It only takes up 210*190 mm of space.
- The microplate is installed vertically into the rotor by the slot on top of the centrifuge. Under the action of surface tension, the liquid is kept on the bottom of the microplate, so we don't need to worry about liquid leakage.
- If you press the button at the front of the centrifuge, the rotor instantly accelerates to 2,500 rpm, but if you release it, the rotor stops.

SPECIFICATIONS

Model	BET-2801
Capacity	2 96-well PCR plates
Power supply	110V-220V/50-60HZ
Power	20W
RPM	2,500r/min
Centrifugal Force	500xg
Rotor	Vertically installed
Standard Configuration	Contains 2 96-well PCR plates
Dimensions (LxWxH)	210mmx190mmx190mm
Weight	3kg
Alt Name	Mini Microplate Centrifuge



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com