



LABORATORY LOW SPEED CENTRIFUGE BET-3401

LABORATORY LOW SPEED CENTRIFUGE BET-3401

OIL CENTRIFUGE

Designed in line with the Test Method for Water and Sediment in Crude Oil by the Centrifuge Method (GB/T6533-2012), this centrifuge is an ideal moisture determination instrument for the oil exploitation industry and R&D institutions, which measures the water and sediment in crude oil by the centrifuge method.



SPECIFICATIONS

Model	BET-3401
Max RCF	2,960xg
Max speed	4,000r/min
Speed accuracy	±30r/min (customizable, multiples of 10)
Motor	High torque variable frequency AC/DC motor, microcomputer control
Temperature control range	
Temperature control accuracy	
Timer range	1min~99h59min59s, with continuous centrifugation and instantaneous centrifugation
Power supply	AC220V50Hz
Noise	≤60dB(A)
Weight	50kg
Dimensions (L x W x H)	450mmx590mmx365mm
Alt Name	Oil Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Description
1	24x15mL	Max RPM: 4,000r/min Max RCF: 2,960xg



1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com