



LABORATORY LOW SPEED CENTRIFUGE BET1BB1 (BCMI-201)

LABORATORY LOW SPEED CENTRIFUGE BET1BB1

HANDHELD CENTRIFUGES



- Delicate in appearance and small in volume, easy to operate and highly stable.
- Patented mute technology allows the machine to operate at ultra-low noise.
- Micro DC motors and carbon fiber rotors.

SPECIFICATIONS

Model	BET1BB1
Old Model	BCMI-201
Max speed	4,000r/min
Max RCF	900xg
Rotor configuration	6x1.5/0.5/0.2mL
Dimensions (LxWxH)	170mmx150mmx130mm
Weight	1.25kg
Alt Name	Handheld Centrifuges



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com