



## LABORATORY HIGH SPEED CENTRIFUGE BET1BB3

# LABORATORY HIGH SPEED CENTRIFUGE BET1BB3

## HANDHELD CENTRIFUGES



- Delicate in appearance and small in volume, easy to operate and highly stable.
- Patented mute technology allows the machine to operate at ultra-low noise.
- Micro DC motors and carbon fiber rotors.

## SPECIFICATIONS

Model	BET1BB3
Max speed	10,000r/min
Max RCF	5,000xg
Rotor configuration	6x1.5/0.5/0.2mL2x8x0.2mL
Dimensions (LxWxH)	170mmx150mmx130mm
Weight	1.25kg
Alt Name	Handheld Centrifuges



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)