



LABORATORY LOW SPEED CENTRIFUGE BET1BB4

LABORATORY LOW SPEED CENTRIFUGE BET1BB4

HANDHELD CENTRIFUGES



- Delicate in appearance and small in volume, easy to operate and highly stable.
- All-steel housing, tailored aluminum alloy rotors, low noise, and small vibrations.
- Timing function.

SPECIFICATIONS

Model	BET1BB4
Max speed	4,000r/min
Max RCF	1,430xg
Rotor configuration	6x20mL
Dimensions (LxWxH)	250(Φ)x270(h)mm
Weight	2.5kg
Alt Name	Handheld Centrifuges



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com