



NATURAL CONVECTION OVEN BEV1W1 (BONC-308)

NATURAL CONVECTION OVEN BEV1W1

NATURAL CONVECTION OVEN

With easy to clean interiors and natural convection circulation system for efficient heating of the chamber, this oven can be relied on for general laboratory applications.



Stainless steel inner chamber without scroll rolled
Silicone rubber door seal
Strong sealing
Nice design

SPECIFICATIONS

Model	BEV1W1
Old Model	BONC-308
Internal dimension/capacity (HxWxD) mm / liter	300x275x280
External dimension (HxWxD) mm	450x635x465
Package dimension (HxWxD) mm	570x720x545
Voltage (V)	220
Power (W)	900
Temp. Range (°C)	50~200
Vacuum degree (Pa)	<133
Temp. Precision (°C)	±0.5
Floating Range (°C)	±1%
Weight (KG) (Wooden shelves)	52.5
Price (USD)	330
Note: Digital display intelligent temperature controller, the timing range is 0~999 minutes	
Alt Name	Vacuum drying oven

APPLICATIONS

Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com