



NATURAL CONVECTION OVEN BEV1W2 (BONC-310)

NATURAL CONVECTION OVEN BEV1W2

NATURAL CONVECTION OVEN

With easy to clean interiors and natural convection circulation system for efficient heating of the chamber, this oven can be relied on for general laboratory applications.



Stainless steel inner chamber without scroll rolled
Silicone rubber door seal
Strong sealing
Nice design

SPECIFICATIONS

Model	BEV1W2
Old Model	BONC-310
Internal dimension/capacity (HxWxD) mm / liter	415x370x345
External dimension (HxWxD) mm	570x750x515
Package dimension (HxWxD) mm	690x830x600
Voltage (V)	220
Power (W)	1000
Temp. Range (°C)	50~200
Vacuum degree (Pa)	<133
Temp. Precision (°C)	±0.5
Floating Range (°C)	±1%
Weight (KG) (Wooden shelves)	75
Price (USD)	460
Note: Digital display intelligent temperature controller, the timing range is 0~999 minutes	
Alt Name	Vacuum drying oven

APPLICATIONS

Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com