



LABORATORY HIGH SPEED CENTRIFUGE BEW1B1 (BBRC-105)

LABORATORY HIGH SPEED CENTRIFUGE BEW1B1

FLOOR-STANDING HIGH-CAPACITY CENTRIFUGE



The max capacity of the swing rotor is 4X1000ml.
It has micro-computer control, electronic lock and AC variable-conversion motor which delivers high torque.
Digital LCD screens make visualized operation convenient.
Imported compressor unit with Non-CFC refrigerant R404a.
It has self-locking, over-speed detection, over-temperature detection, imbalance detection that is an automatic alarm device.
Essential configuration included rotor exchange, speed/RCF, time, temperature, 9 levels acceleration and deceleration is the most economical choice.
Silastic sealed ring comfort to GMP, US FDA certificated.

SPECIFICATIONS

Model	BEW1B1
Old Model	BBRC-105
Max speed	5000 r/min
Max RCF	4100 xg
Max capacity	4x1000 ml
Speed accuracy	±10 r/min
Time setting	1s~99h 59min59s
Temperature setting	-20°C ~ +40°C
Temperature accuracy	±1°C
Compressor unit	Yes
Noise	<65 dB(A)
Power supply	AC220V, 50Hz-60Hz
Consumption	2 kW
Alt Name	Floor-standing High-Capacity Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Rotor	Capacity: 6x300ml Max Speed: 5000r/min Max RCF: 3850xg Adaptor: 50ml, 100ml
2	Rotor	Capacity: 4x500ml Max Speed: 4000r/min Max RCF: 3780xg
3	Adapter	500 ml Bucket
4	Adapter	4x20x1.5ml
5	Adapter	4x20x5ml / 7ml Blood Vacuum Tube
6	Adapter	4x15x10ml Internal Thread
7	Adapter	4x8x15ml
8	Adapter	500 ml PE Bottle
9	Adapter	4x9x20ml
10	Adapter	4x4x50ml

11	Adapter	4x2x100ml
12	Adapter	4x250ml
13	Rotor	Capacity: 4x750ml Max Speed: 4000r/min Max RCF: 3780xg
14	Adapter	750ml bucket
15	Adapter	750 ml PE Bottle
16	Rotor	Capacity: 4x1000ml Max Speed: 4000r/min Max RCF: 4100xg
17	Adapter	1000ml bucket
18	Adapter	4x21x1.5ml
19	Adapter	4x24x5ml / 7ml Blood Vacuum Tube
20	Adapter	4x37x5ml / 7ml G RIA Tube
21	Adapter	4x20x10ml Internal Thread
22	Adapter	4x14x15ml
23	Adapter	1000 ml PE Bottle
24	Adapter	4x12x20ml
25	Adapter	4x5x50ml
26	Adapter	4x3x100ml
27	Adapter	4x250ml
28	Adapter	4x500ml
29	Rotor	Capacity: 6x1000ml Max Speed: 8000r/min Max RCF: 15190xg Adaptor: 50ml, 100ml, 250ml





25



26



27



28



29



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com