



SHAKING WATER BATH BEX1D1 (BBSW-101)

SHAKING WATER BATH BEX1D1

SHAKING WATER BATH



Microcomputer temperature and frequency controller, LCD display screen;
Stainless steel chamber;
With power outage recovery and parameter memory function;
Universal spring fixture;
Equipped with timing function.

SPECIFICATIONS

Model	BEX1D1
Old Model	BBSW-101
Type	Water Bath
Temp. Range [°C]	RT+5~100
Oscillation Mode	Orbital
Rotation Rate [rpm]	0~220
Vibration Amplitude [mm]	20
Power [W]	1600
Timer Range [min/h]	0~9999
Internal Dimensions WxDxH [mm]	480x390x170
Housing Dimensions WxDxH [mm]	705x515x470
Tray Size WxDxH [mm]	380x300x170
Full Load Configuration [mlxpcs]	500x6 or 1000x3
Net Weight/Gross weight [kg]	45/61
Package Size [mm]	810x620x600
Alt Name	Shaking Water Bath



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com