



## MINI CENTRIFUGE BFA-1805

# MINI CENTRIFUGE BFA-1805

## MINI CENTRIFUGE

Mini Centrifuge is perfectly designed with novel and unique appearance. It is pretty and skillful, equipped with 2 types of rotors and several types of tube carriers for 2.0ml / 0.5ml / 0.2ml tubes and 0.2ml x 8 PCR tube strips. It is designed user-friendly. Open/close the cover will automatically start/stop the centrifuge. Mini centrifuge series are built-in timing function.



Large radius rotor with 8 tube holes, centrifugal force is 1.7times bigger than ordinary 6-hole rotor at the same rolling speed.

The cover press-button is outward designed with start/stop function. Press the button, the cover will open automatically to 95°, which is very convenient for one-hand operation.

The rotor is clip-on fixed. Very convenient for rotors exchanging yet fixed solid reliable.

LED displays speed and timing.

Working is silent and stable. Free-maintenance motor, long-life and safe.

Smart appearance. Multi-purpose to meet different experimental demands.

## SPECIFICATIONS

|                |                                                                                                                   |
|----------------|-------------------------------------------------------------------------------------------------------------------|
| Model          | BFA-1805                                                                                                          |
| Speed Range    | 3krpm~10krpm (increment:1000rpm)                                                                                  |
| RCF            | approx.400g~4700g                                                                                                 |
| Rotor capacity | 8x1.5/2.0ml, 8x0.5ml, 2x8x0.2ml strips(With 0.5ml and 0.2ml tube adapter)2x8x0.2ml PCR Strips(when speed≤6000rpm) |
| Time Range     | 1s ~9999min or continuous                                                                                         |
| Working Noise  | ≤55dB                                                                                                             |
| Voltage        | AC 220V / AC 110V, 50/60Hz                                                                                        |
| Power          | ≤45W                                                                                                              |
| Fuse           | 250V, 1A, Φ5x20                                                                                                   |
| Dimension      | W.156 x D.176 x H.121mm                                                                                           |
| Net Weight     | 1.5kgs                                                                                                            |
| Alt Name       | Mini Centrifuge                                                                                                   |

## ITEMS INCLUDED

| No | Name           |
|----|----------------|
| 1  | Combined rotor |
| 2  | Tube adapters  |



1



2



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)