



MULTISTATION STIRRER BFA-6401

MULTISTATION STIRRER BFA-6401

MULTI-POINT ULTRA-THIN MAGNETIC STIRRER

The multi-thin ultrathin magnetic stirrer can precisely and stably stir the liquid in a wide range, especially suitable for small volume samples, and is an essential tool for petroleum, chemical, environmental protection, education and scientific research.



Ultra-light and ultra-thin, novel appearance, compact size, saving experiment space.

The controller and the agitator are separated. The number and combination of stirrer can be adjusted arbitrarily, up to 6 stirrer.

Adopt the connection method of USB interface, easy to plug and unplug. Each stirrer has its own speed and time, which makes it easier to compare experiments.

Touch keys and sliding touch adjustment, durable, novel and fashionable.

Advanced non-wearing magnetic coil technology, the operation is quieter and more stable, and the service life is longer.

The LCD screen displays the parameters in real time, the LED indicator shows the working status, and the instrument running status is more intuitive.

SPECIFICATIONS

Model	BFA-6401
Type	Controller
Max. stirring volume (H2O)	-
Motor output power	-
Speed setting (running) range	20-1500rpm (increment 10)
Timing range	0-9 hours (0 means infinity, increment 5min)
Display method	LCD
Working plate material	-
Working plate size	-
Dimensions	W.165xD.105xH.45mm
Weight	0.48kgs
Power adapter input	AC100-240V 50/60Hz 60VA
Power adapter output	DC 24V, 2A
Alt Name	Multi-point ultra-thin magnetic stirrer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com