



MULTISTATION STIRRER BFA-6402

MULTISTATION STIRRER BFA-6402

MULTI-POINT ULTRA-THIN MAGNETIC STIRRER

The multi-thin ultrathin magnetic stirrer can precisely and stably stir the liquid in a wide range, especially suitable for small volume samples, and is an essential tool for petroleum, chemical, environmental protection, education and scientific research.



Ultra-light and ultra-thin, novel appearance, compact size, saving experiment space.

The controller and the agitator are separated. The number and combination of stirrer can be adjusted arbitrarily, up to 6 stirrer.

Adopt the connection method of USB interface, easy to plug and unplug. Each stirrer has its own speed and time, which makes it easier to compare experiments.

Touch keys and sliding touch adjustment, durable, novel and fashionable. Advanced non-wearing magnetic coil technology, the operation is quieter and more stable, and the service life is longer.

The LCD screen displays the parameters in real time, the LED indicator shows the working status, and the instrument running status is more intuitive.

SPECIFICATIONS

| | |
|-------------------------------|---|
| Model | BFA-6402 |
| Type | Stirrer |
| Max. stirring volume (H2O) | 1L |
| Motor output power | 4W |
| Speed setting (running) range | (20-1500rpm) |
| Timing range | - |
| Display method | - |
| Working plate material | PVC |
| Working plate size | Φ120mm |
| Dimensions | W.120xD.120xH.19mm |
| Weight | 0.18kgs |
| Power adapter input | - |
| Power adapter output | - |
| Alt Name | Multi-point ultra-thin magnetic stirrer |



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com