



ORBITAL SHAKER BFA1BA1 (BSOT-205)

ORBITAL SHAKER BFA1BA1

DECOLORIZATION SHAKER

Decolorization shaker is widely used in coomassie blue dyeing and decolorization shaking, electrophoresis gel fixing, silver nitrate dyeing fixing, dyeing, development and so on.



Compact appearance, space saving, powerful function and strong compatibility. Suction cup type machine feet, super shockproof, high speed and smooth, no noise.

Large platform design, placing more samples, strong load capacity.

Three models are available, the configuration performance is from low to high, to meet more choices.

SPECIFICATIONS

Model	BFA1BA1
Old Model	BSOT-205
Speed Range	0~300rpm
Time Range	/
Digital Display	/
Radius of gyration	10mm
Max Load Capacity	3Kgs
Motor Parameters	DC motor stepless speed regulation
Platform Size	330x430mm
Beep Alarm	No
Voltage	AC110/220V,50/60Hz
Power	25W
Fuse	250V,1A,φ5x20
Dimension	W:330xD:490xH:145mm
Net Weight	9.5kgs
Alt Name	Decolorization shaker

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Platform TS-A	Stainless steel platform and rubber mat combination
2	Conical Flask Stand TS-B	250ml can put 12pcs
3	Conical Flask Stand TS-C	100ml can put 20pcs
4	Universal Platform TS-D	



1



2



3



4



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com