



## BENCHTOP SHAKING INCUBATOR BFA1BZ1 (BDIB-112)

# BENCHTOP SHAKING INCUBATOR BFA1BZ1

## MICRO-PLATE SHAKER INCUBATOR

Micro-plate incubator is a well-designed micro processing techniques incubator which well combined with PID controlling. It is mainly used for solution incubation of ELISA plate (96-well/384-well), cell culture plate (24-well / 48-well / 96-well) at proper temperature.



LCD display. Setting temperature and actual temperature simultaneously displayed. Diminishing time displayed.

Friendly human-machine operation interface. Streamline shape, pretty and elegant. Convenient for cleaning.

Uniformly upper & down heating for samples in micro plates. The experiment effect is exact correspondence of the setting and actual temperature.

Microprocessors controls temperature and time, good temperature control linear, small fluctuations.

2/4 standard ELISA plates and microplates can be placed, and sound alarm signal after the program finished.

## SPECIFICATIONS

Model	BFA1BZ1
Old Model	BDIB-112
Temp. Control Range	R.T.+5°C ~70°C
Temp. Setting Range	5°C~70°C
Temp. Stability@37°C	±0.2°C
Platform Temp. Uniformity	±0.5°C
Temp. Display Accuracy	0.1°C
Heating Speed	<25min (20°C to 70°C)
Time Range	1min~99h59min
Sample Capacity	2 Pieces of Plates(Height <25mm)
Voltage	AC 220V/AC 110V, 50/60Hz
Power	150W
Fuse	250V, 1A/2A, Φ5x20
Dimension	W.284 x D.264 x H.157mm
Net Weight	3.3kgs
Alt Name	Micro-plate Shaker Incubator



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)