



HORIZONTAL GEL ELECTROPHORESIS BFA1CT3 (BGEL-401)

HORIZONTAL GEL ELECTROPHORESIS BFA1CT3

HORIZONTAL GEL ELECTROPHORESIS

Horizontal Gel Electrophoresis is a one-body electrophoresis Unit, which has compact design, stable performance, easy operation and good reproducibility and cost performance. It is widely used in chemical field, ideal for teaching, scientific research, food and medical treatment. Horizontal gel electrophoresis unit is mainly used for separation and test of nucleic acid.



New and unique appearance, mini and portable; The electrophoresis tank and power can be separated.

The cover with safety design, flip cover immediately power off, to ensure safety; Transparent cover and open hole design, easy to heat dissipation and observation;

Top handle design, easy to pick and place.

With seven output voltages, it can be called by one button, which is convenient and fast, and saves standby consumption.


Specially designed Electrophoresis shape, reducing the temperature difference between the buffer, while the heat benefit is improved. Put the gel into a limited position without drift. With the pause function, timing function, can leave after press the start button, with Alarm function, buzzer alert at the end.

SPECIFICATIONS

Model	BFA1CT3
Old Model	BGEL-401
Capacity	Max.500ml
Inner diameter dimension	W.150xD.150xH.47mm
Multi-channel suitability	Suitable for multi-channel pipettor
Timing function	1min~99h59min or contin., with pause function
Security Settings	The power is automatically paused when the cover is lifted
Storage function	Automatically memorize the previous operating voltage and time
Input power	AC110V, 50/60Hz
Output voltage	Constant peak voltage (140V) , impulse control
Dimensions	W.260xD.170xH.68mm
Optional output voltage	135V 100V 70V 50V 35V 25V 18V
Net weight	0.78kgs
Alt Name	Horizontal Gel Electrophoresis

ACCESSORIES FOR PURCHASE

No	Name	Description	Specifications/dimensions	Quantity
1		Small size plate for gel	60mm(W)x60mm(L)	2
2		Medium size plate for gel	120mm(W)x60mm(L)	1
3		big size plate for gel	120mm(W)x120mm(L)	1
4	Gel Comb 1	Gel Comb 1 1.0mm thickness 25/11 teeth (Can used multichannel gun to add sample)		4
5	Gel Comb 2	Gel Comb 2 1.5mm thickness 18/8 teeth		1
6	Gel Comb 3	Gel Comb 3 1.5mm thickness 13/6 teeth		1
7	Gel Comb 4	Gel Comb 4 2.0mm thickness 3/2 teeth		1
8		Base for gel making (double side)	138mm(W)x138mm(L)x36mm(H)	1



FEATURES

New and unique appearance, mini and portable; The electrophoresis tank and power can be separated.

The cover with safety design, flip cover immediately power off, to ensure safety; Transparent cover and open hole design, easy to heat dissipation and observation;

Top handle design, easy to pick and place.

With seven output voltages, it can be called by one button, which is convenient and fast, and saves standby consumption.

Specially designed Electrophoresis shape, reducing the temperature difference between the buffer, while the heat benefit is improved. Put the gel into a limited position without drift. With the pause function, timing function, can leave after press the start button, with Alarm function, buzzer alert at the end.

Highlight digital tube displays the set parameters, default remember the output voltage and time Settings of the previous operation.

When the electrophoresis current exceeds the peak current, will automatically adjust the voltage to reduce the current value, improve the experimental reproducibility.

Equipped with variety of specifications of gel trays / combs, compatible with a variety of experimental needs, compatible with high-throughput electrophoresis and multi-channel pipette operation.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com