



MINI CENTRIFUGE BFA1S5

MINI CENTRIFUGE BFA1S5

MINI CENTRIFUGE

Mini Centrifuge is perfectly designed with novel and unique appearance. It is pretty and skillful, equipped with 2 types of rotors and several types of tube carriers for 2.0ml / 0.5ml / 0.2ml tubes and 0.2ml x 8 PCR tube strips. It is designed user-friendly. Open/close the cover will automatically start/stop the centrifuge. Mini centrifuge series are built-in timing function.



Large radius rotor with 8 tube holes, centrifugal force is 1.7times bigger than ordinary 6-hole rotor at the same rolling speed.

The cover press-button is outward designed with start/stop function. Press the button, the cover will open automatically to 95°, which is very convenient for one-hand operation.

The rotor is clip-on fixed. Very convenient for rotors exchanging yet fixed solid reliable.

LED displays speed and timing.

Working is silent and stable. Free-maintenance motor, long-life and safe.

Smart appearance. Multi-purpose to meet different experimental demands.

SPECIFICATIONS

Model	BFA1S5
Speed Range	3krpm~10krpm (increment:1000rpm)
RCF	approx.400g~4700g
Rotor capacity	8x1.5/2.0ml, 8x0.5ml, 2x8x0.2ml strips(With 0.5ml and 0.2ml tube adapter)2x8x0.2ml PCR Strips(when speed≤6000rpm)
Time Range	1s ~9999min or continuous
Working Noise	≤55dB
Voltage	AC 220V / AC 110V, 50/60Hz
Power	≤45W
Fuse	250V, 1A, Φ5x20
Dimension	W.156 x D.176 x H.121mm
Net Weight	1.5kgs
Alt Name	Mini Centrifuge

ACCESSORIES FOR PURCHASE

No	Name
1	Combined rotor
2	Tube adapters



1



2



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com