



## FLOOR TYPE FREEZE DRYER BFB1BD1 (BFFT-315-A)

# FLOOR TYPE FREEZE DRYER BFB1BD1

## FREEZE DRYER

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.



## SPECIFICATIONS

|                       |                                      |
|-----------------------|--------------------------------------|
| Model                 | BFB1BD1                              |
| Old Model             | BFFT-315-A                           |
| Condenser Temperature | -110°C                               |
| Freeze Drying area    | 0.18 or 0.27 m <sup>2</sup>          |
| Bulk capacity         | 1.8 or 2.7 L                         |
| Tray                  | φ240 mm, 4 or 6 pcs                  |
| Alt Name              | Standard type electric heating shelf |

## APPLICATIONS

Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)