



FLOOR TYPE FREEZE DRYER BFB1BF1 (BFFT-317-A)

FLOOR TYPE FREEZE DRYER BFB1BF1

FREEZE DRYER

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.



SPECIFICATIONS

Model	BFB1BF1
Old Model	BFFT-317-A
Condenser Temperature	-110°C
Freeze Drying area	0.18 or 0.27 m ²
Bulk capacity	1.8 or 2.7 L
Tray	φ240mm, 4 or 6 pcs
Round bottom flask	1000 ml, 500 ml, 250 ml, 100 ml, 50 ml
Alt Name	Manifold type Electric heating shelf

APPLICATIONS

Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com