



## FLOOR TYPE FREEZE DRYER BFB1BG1 (BFFT-318-A)

# FLOOR TYPE FREEZE DRYER BFB1BG1

## FREEZE DRYER

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.



## SPECIFICATIONS

Model	BFB1BG1
Old Model	BFFT-318-A
Condenser Temperature	-110 °C
Freeze Drying area	0.09 m <sup>2</sup>
Bulk capacity	0.9 L
Tray	φ200 mm, 3 pcs
Round bottom flask	1000 ml, 500 ml, 250 ml, 100 ml, 50 ml
Alt Name	Top Press with Manifold type electric heating shelf

## APPLICATIONS

Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)