



FLOOR TYPE FREEZE DRYER BFB1S2 (BFFT-303-B)

FLOOR TYPE FREEZE DRYER BFB1S2

FREEZE DRYER

Compact, stainless steel unit with casters for easy movement, Floor type freeze dryer provides unmatched process accuracy and reliability in performance. It optimizes your space requirements and reduces energy consumption.



SPECIFICATIONS

Model	BFB1S2
Old Model	BFFT-303-B
Condenser Temperature	-110°C
Freeze Drying area	0.18 m ²
Bulk capacity	1.8 L
Tray	φ200mm, 6pcs
Round bottom flask	1000ml, 500ml, 250ml, 100ml, 50ml
Alt Name	Manifold Type

APPLICATIONS

Perishable Substances, Waste Products, Laboratory, Research, Proteins, Organic Tissues, Waste Products, Plant Material, Polymers, Pharmaceuticals, Nutraceuticals, Plant material.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com