



BULLET ICE MAKER BFC1C1 (BIBU-101)

BULLET ICE MAKER BFC1C1

VERTICAL ICE MAKER



Luxurious appearance and compact design;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Cylindrical bullet ice with elegant appearance.
Ice thickness can be adjusted, with a wider scope of application.
Bullet Ice

SPECIFICATIONS

Model	BFC1C1
Old Model	BIBU-101
Ice making capacity	20kg/24h
Power/voltage	220w
Ice storage capacity	6kg
Ice shape	Bullet shape
N.W./kgs	28kg
G.W./kgs	31kg
Product dimension	330x470x605mm
Packing Dimension	400x550x680
20GP	168
40GP	348
40HQ	464
Alt Name	Vertical Ice Maker

FEATURES

Luxurious appearance and compact design;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Cylindrical bullet ice with elegant appearance.
Ice thickness can be adjusted, with a wider scope of application.



Bullet Ice



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com