



## BULLET ICE MAKER BFC1C4 (BIBU-109)

# BULLET ICE MAKER BFC1C4

## VERTICAL ICE MAKER



Luxurious appearance and compact design;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Cylindrical bullet ice with elegant appearance.  
Ice thickness can be adjusted, with a wider scope of application.  
Bullet Ice

## SPECIFICATIONS

Model	BFC1C4
Old Model	BIBU-109
Ice making capacity	80kg/24h
Power/voltage	580w
Ice storage capacity	25kg
Ice shape	Bullet shape
N.W./kgs	51kg
G.W./kgs	56kg
Product dimension	500x612x900mm
Packing Dimension	580x700x970mm
20GP	64
40GP	136
40HQ	136
Alt Name	Vertical Ice Maker

## FEATURES

Luxurious appearance and compact design;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Cylindrical bullet ice with elegant appearance.  
Ice thickness can be adjusted, with a wider scope of application.



Bullet Ice



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)