



BULLET ICE MAKER BFC1C5 (BIFL-212)

BULLET ICE MAKER BFC1C5

VERTICAL ICE MAKER



Luxurious appearance and compact design;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Cylindrical bullet ice with elegant appearance.
Ice thickness can be adjusted, with a wider scope of application.
Bullet Ice

SPECIFICATIONS

Model	BFC1C5
Old Model	BIFL-212
Ice making capacity	120kg/24h
Power/voltage	760w
Ice storage capacity	50kg
Ice shape	Bullet shape
N.W./kgs	70kg
G.W./kgs	75kg
Product dimension	680x610x1050mm
Packing Dimension	750x680x1100mm
20GP	48
40GP	96
40HQ	96
Alt Name	Vertical Ice Maker

FEATURES

Luxurious appearance and compact design;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Cylindrical bullet ice with elegant appearance.
Ice thickness can be adjusted, with a wider scope of application.



Bullet Ice



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com