



## CUBE ICE MAKER BFC1D2 (BIFL-214)

# CUBE ICE MAKER BFC1D2

## SPLIT-TYPE ICE MAKER



Split-type structure, large ice bin, large ice storage capacity and strong applicability;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Square ice cubes with beautiful shapes.  
Square

## SPECIFICATIONS

Model	BFC1D2
Old Model	BIFL-214
Ice Making Capacity	300kg/24h
Ice Storage Capacity	150kg
Ice Shape	Square
Power/Voltage (W/V)	1240w/220v
Cooling Type	Water cooling
Product Dimension (mm)	760x850x1620
Machine Appearance	Split machine
Alt Name	Split-Type Ice Maker

## FEATURES

Split-type structure, large ice bin, large ice storage capacity and strong applicability;  
Fully automatically computer controlled, convenient for use;  
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;  
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;  
Square ice cubes with beautiful shapes.



Square



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)