



BENCHTOP LOW SPEED CENTRIFUGE BFD-10501

LOW SPEED CENTRIFUGE BFD-10501

BENCHTOP LOW SPEED CENTRIFUGE



SPECIFICATIONS

Model	BFD-10501
Max. Speed	4200rpm
Max. RCF	2680xg
Max. Capacity	12 x 15 ml
Running Time	1 min-99 min
Drive System	Direct, induction motor drive (brushless)
Speed Accuracy	±20 r/min
Accel/Decel Profiles	9/9
Noise Level	≤54 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	200 W
Dimension (W x D x H)	13.4 x 17.3 x 11.4 in / 34 x 44 x 29 cm
Weight w/o Rotor	40 lb / 18 kg
Alt Name	Benchtop Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Fixed-Angle Rotor	8 x 15 ml (conical)	4 200 rpm	2 190 Xg
2	adapters			
3	adapters			
4	adapters			
5	adapters			
6	adapters			
7	adapters			
8	No. 2 Fixed-Angle Rotor	12 x 15 ml (round)	4 200 rpm	2 680 Xg
9	adapters			
10	adapters			
11	No. 3 Fixed-Angle Rotor	4 x 50 ml		
12	adapters			
13	adapters			
14	adapters			
15	adapters			

16	adapters			
17	adapters			
18	adapters			
19	adapters			



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com