



BENCHTOP LOW SPEED CENTRIFUGE BFD-10601

LOW SPEED CENTRIFUGE BFD-10601

BENCHTOP LOW SPEED CENTRIFUGE



SPECIFICATIONS

Model	BFD-10601
Max. Speed	4950rpm
Max. RCF	3050xg
Max. Capacity	24 x 10 ml blood collection/urine tube
Running Time	1 min-99 min
Drive System	Direct, induction motor drive (brushless)
Speed Accuracy	±20 r/min
Accel/Decel Profiles	9/9
Noise Level	≤60 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	200 W
Dimension (W x D x H)	13.4 x 17.3 x 11.4 in / 34 x 44 x 29 cm
Weight w/o Rotor	40 lb / 18 kg
Alt Name	Benchtop Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Fixed-Angle Rotor	8 x 15 ml (conical)	4 950 rpm	3 050 Xg
2	adapters			
3	adapters			
4	adapters			
5	adapters			
6	adapters			
7	adapters			
8	No. 2 Fixed-Angle Rotor	24 x 10 ml blood collection/urine tube	4 200 rpm	2 680 Xg
9	adapters			
10	adapters			
11	No. 3 Fixed-Angle Rotor	12 x 15 ml (round)		
12	adapters			
13	adapters			
14	No. 4 Fixed-Angle Rotor	4 x 50 ml		

15	adapters			
16	adapters			
17	adapters			
18	adapters			
19	adapters			
20	adapters			
21	adapters			
22	adapters			



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com