



BENCHTOP LOW SPEED CENTRIFUGE BFD-10701

LOW SPEED CENTRIFUGE BFD-10701

BENCHTOP LOW SPEED CENTRIFUGE



SPECIFICATIONS

Model	BFD-10701
Max. Speed	4200rpm
Max. RCF	2680 xg
Max. Capacity	24 x 10 ml blood collection/urine tube
Running Time	1 min-99 min
Drive System	Direct, induction motor drive (brushless)
Speed Accuracy	±20 r/min
Accel/Decel Profiles	9/9
Noise Level	≤60 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	250 W
Dimension (W x D x H)	11.8 x 14.8 x 10.2 in / 30 x 37.5 x 26 cm
Weight with Rotor	40 lb / 18 kg
Alt Name	Benchtop Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Fixed-Angle Rotor	24 x 10 ml blood collection/urine tube	4 200 rpm	2 680 Xg
2	adapters			
3	adapters			



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com