



BENCHTOP LOW SPEED CENTRIFUGE BFD-10801

LOW SPEED CENTRIFUGE BFD-10801

BENCHTOP LOW SPEED CENTRIFUGE



SPECIFICATIONS

Model	BFD-10801
Max. Speed	5500rpm
Max. RCF	4300xg
Max. Capacity	12 x 15 ml
The shortest time of acceleration to the max speed	≤10 s
Control System	Microprocessor
Drive System	Direct, induction motor drive
Display	Full color touchscreen interface
Programmability	Via easy access touchscreen menu
Multi-stage Centrifugal Mode	Max up to 5 stages of speed & time
User Lockout Function	Yes
Data Communication	USB port (Run history can be output through USB port)
Program Storage	1,000 user-defined
Running Time	1 s-99 h:59 min:59 s, Short
Timer	Count up or down
Speed Accuracy	±10 r/min
Accel/Decel Profiles	40/40
Noise Level	≤60 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	500 W
Dimension (W x D x H)	15.7 x 19.3 x 12.2 in / 40 x 49 x 31 cm
Weight with Rotor	60 lb / 27 kg
Alt Name	Benchtop Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Fixed-Angle Rotor	12 x 15 ml	5 500 rpm	4 300 Xg
2	adapters			
3	adapters			
4	adapters			
5	adapters			

6	adapters			
7	adapters			



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com