



BENCHTOP LOW SPEED CENTRIFUGE BFD-10901

LOW SPEED CENTRIFUGE BFD-10901

BENCHTOP LOW SPEED CENTRIFUGE



SPECIFICATIONS

Model	BFD-10901
Max. Speed	5000rpm
Max. RCF	4650xg
Max. Capacity	8 x 100 ml
Running Time	1 min - 99 min 59 s, Pulse
Drive System	Direct, induction motor drive (brushless)
Speed Accuracy	±10 r/min
Programs	100 user-defined
Accel/Decel Profiles	20/20
Noise Level	≤60 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	500 W
Dimension (W x D x H)	17 x 21.7 x 13.8 in / 43 x 55 x 35 cm
Weight w/o Rotor	88 lb / 40 kg
Alt Name	Benchtop Low Speed Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Swing-out Rotor		5 000 rpm	4 650 Xg
2	Rotor Cross			
3	Tube Racks	4 x 50 ml		
4	Tube Racks	4 x 100 ml		
5	No. 2-1 Swing-out Rotor		4 000 rpm	2 980 Xg
6	Rotor Cross			
7	Tube Racks	8 x 50 ml		
8	Tube Racks	8 x 100 ml		
9	Tube Racks	24 x 10 ml		
10	Tube Racks	24 x 15 ml		
11	adapters			
12	adapters			
13	adapters			
14	adapters			

15	adapters			
16	adapters			
17	Tube Racks	32 x 10 ml		
18	Tube Racks	32 x 15 ml		
19	adapters			
20	adapters			
21	adapters			
22	adapters			
23	adapters			
24	adapters			
25	No. 2-2 Swing-out Rotor	48 x 1.6 -7 ml blood collection / urine tube		
26	No. 3 Swing-Bucket Rotor	2 x 3 pcs microplates (96 wells) or 2 x 3 pcs microtiter plates or 2 x 2 pcs cell culture plates or 2 x 1 pc deep well plates	4 000 rpm	1 970 Xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com