



MINI CENTRIFUGE BFD-11501

MINI CENTRIFUGE BFD-11501

MINI CENTRIFUGE



SPECIFICATIONS

Model	BFD-11501
Max. Speed	3000rpm
Max. RCF	600xg
Max. Capacity	2 x 96 well PCR plates
Running Time	1 s-99 s, 1 min-99 min, Pulse
Drive System	DC 24 V, Brushless DC motor
Power Supply (power adapter)	AC 100-240 V, 50/60 Hz
Dimension (W x D x H)	7.9 x 9 x 7.5 in / 20 x 23 x 19 cm
Weight	5.7 lb / 2.6 kg
Alt Name	Mini Centrifuge

ITEMS INCLUDED

No	Name	Capacity	Speed	RCF
1	Rotor	2 x 96 well PCR plates	3000rpm	600xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com