



MICRO CENTRIFUGE BFD-11701

MICRO CENTRIFUGE BFD-11701

MINI CENTRIFUGE



SPECIFICATIONS

Model	BFD-11701
Max. Speed	10000rpm
Max. RCF	6100xg
Max. Capacity	
Rotor Capacity	8 x 1.5/2.0 ml + 6 x 0.5 ml + 2 x PCR-Strips / 16 x 0.2 ml (Three capacities in one rotor)
Running Time	1 s - 99 min : 59 s
Drive System	Brushless DC motor
More Features	Lid open protection; Auto power motor off when lid is open; Stop running within 3 seconds when lid is open
Noise Level	≤57 dB(A)
Power Supply	AC 100-240 V, 50/60 Hz; Lithium battery pack optional for outdoor centrifugation
Dimension (W x D x H)	6.7 x 9.9 x 5.1 in / 17 x 20 x 13 cm
Weight without Lithium Battery Pack	2.6 lb / 1.2 kg
Alt Name	Mini Centrifuge

ITEMS INCLUDED

No	Name	Capacity	Speed	RCF
1	Rotor	8 x 1.5/2.0 ml + 6 x 0.5 ml + 2 x PCR-Strips / 16 x 0.2 ml	10000 rpm	6100 xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com