



MICRO CENTRIFUGE BFD-11801

MICRO CENTRIFUGE BFD-11801

MINI CENTRIFUGE



SPECIFICATIONS

Model	BFD-11801
Max. Speed	12000rpm
Max. RCF	9500xg
Max. Capacity	
Rotor Capacity	12 x 1.5/2.0/0.5 ml + 4 x PCR-Strips / 32 x 0.2 ml (Two capacities in one rotor)
Alternate Rotor Capacity	4 x 12 x 0.2 ml / 4 x PCR-Strips
Running Time	1 s - 99 min : 59 s
Drive System	DC 24 V, Brushless DC motor
Power Supply	AC 100-250 V, 50/60 Hz
More Features	Max. accelerate from 0 rpm to max speed within 12 seconds; Two stop programs: 1) Press key "Open" for unbraked deceleration (≤ 15 s); 2) Stop running within 3 seconds when lid is open
Noise Level	≤ 42 dB(A)
Dimension (W x D x H)	7.9 x 9 x 4.7 in / 20 x 23 x 12 cm
Weight with Rotor	4.4 lb / 2 kg
Alt Name	Mini Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Fixed-Angle Rotor	12 x 1.5/2.0/0.2/0.5 ml + 4 x PCR-Strips	12 000 rpm	9 500xg
2	No. 2 Fixed-Angle Rotor	4 x 12 x 0.2ml or 4 x PCR-Strips	12 000 rpm	5 960 Xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com