



MICRO CENTRIFUGE BFD-11901

MICRO CENTRIFUGE BFD-11901

MINI CENTRIFUGE



SPECIFICATIONS

Model	BFD-11901
Max. Speed	15000 rpm
Max. RCF	15093xg
Max. Capacity	8 x 5 ml
Running Time	1 s - 99 min : 59 s, Short
Drive System	Direct, brushless motor drive
Speed Accuracy	±10 r/min
Accel/Decel Profiles	3/3
More Features	Max. accelerate from 0 rpm to 15,000 rpm within 10 seconds; Max. accelerate from 15,000 rpm to 0 rpm within 10 seconds
Max. Power Consumption	100 W
Noise Level	≤56 dB(A)
Power Supply	AC 100-250 V, 50/60 Hz, 3A
Dimension (W x D x H)	8.7 x 9.5 x 6.7 in / 22 x 24 x 17 cm
Weight without Lithium Battery Pack	8.4 lb / 3.8 kg
Alt Name	Mini Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Fixed-Angle Rotor	12 x 1.5/2.0 ml (with adapters for 0.2/0.5ml tube)	15 000 rpm	15 093 Xg
2	No. 2 Fixed-Angle Rotor	8 x 5 ml	15 000 rpm	15 093 Xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com