



PRP CENTRIFUGE BFD-12201

PRP CENTRIFUGE BFD-12201

PRF CENTRIFUGE



SPECIFICATIONS

Model	BFD-12201
Max. Speed	3500rpm
Max. RCF	1640xg
Max. Capacity	12 x 10 ml blood collection tube
Drive System	Direct induction motor drive
Speed Accuracy	± 20 r/min
Noise Level	≤ 55 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	200 W
Dimension (W x D x H)	13.4 x 17 x 13 in / 34 x 43 x 33 cm
Weight with Rotor	40 lb / 18 kg
Alt Name	PRP Centrifuge

ITEMS INCLUDED

No	Name	Capacity	Speed	RCF
1	Rotor	12 x 10 ml blood collection tube	3500rpm	1640xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com