



CYTO CENTRIFUGE BFD-12301

CYTO CENTRIFUGE BFD-12301

CYTO CENTRIFUGE



SPECIFICATIONS

Model	BFD-12301
Max. Speed	4000rpm
Max. RCF	2170xg
Max. Capacity	10 x 10 ml Cyto chamber
Running Time	1 min-99 min
Drive System	Direct induction motor drive
Speed Accuracy	± 20 r/min
Accel/Decel Profiles	9/9
Noise Level	≤ 58 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	200 W
Dimension (W x D x H)	13.4 x 17 x 13 in / 34 x 43 x 33 cm
Weight w/o Rotor	46 lb / 21 kg
Alt Name	Cyto Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Rotor	10 x 10 ml Cyto chamber	4000rpm	2170xg
2	No. 2 Rotor	6 x 9 ml Cyto chambe	4000rpm	2170xg
3	Consumables for No. 1 Rotor, including Disposable Cyto Chamber. Filter Paper.		4000rpm	2170xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com