



DAIRY TESTING CENTRIFUGE BFD-12401

MILK FAT CENTRIFUGE BFD-12401

DAIRY TESTING CENTRIFUGE



SPECIFICATIONS

Model	BFD-12401
Max. Speed	3300rpm
Max. RCF	1920xg
Max. Capacity	4 x 50 ml Swing-out Rotor for BABCOCK TEST BOTTLES (Tube Size: $\Phi 38 \times 170$ mm)
Running Time	1 min - 9999 min 59 s, Pulse, Continuous
Drive System	Direct, induction motor drive
Speed Accuracy	± 10 r/min
Programs	12 user-defined
Accel/Decel Profiles	19/19
Noise Level	≤ 60 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	500 W
Dimension (W x D x H)	17.7 x 21.7 x 14.4 in / 45 x 55 x 36.5 cm
Weight with Rotor	95 lb / 43 kg
Alt Name	Dairy Testing Centrifuge

ITEMS INCLUDED

No	Name	Capacity	Speed	RCF
1	Swing-out Rotor	4 x 50 ml	3300rpm	1920xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com