



DAIRY TESTING CENTRIFUGE BFD-12601

MILK FAT CENTRIFUGE BFD-12601

DAIRY TESTING CENTRIFUGE



SPECIFICATIONS

Model	BFD-12601
Max. Speed	3300rpm
Max. RCF	1920xg
Max. Capacity	8 x 30 mL Angle Rotor (for Butyrometer)
Running Time	1 min - 99 h : 59 min, Pulse
Drive System	Direct, induction motor drive (Inverter Control)
Speed Accuracy	±10 r/min
Programs	12 user-defined
Accel/Decel Profiles	19/19
Noise Level	≤60 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	1.1 kW
Dimension (W x D x H)	21.7 x 25.6 x 33 in / 55 x 65 x 84 cm
Weight w/o Rotor	293 lb / 133 kg
Alt Name	Dairy Testing Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Speed	RCF
1	No. 1 Fixed-Angle Rotor	8 x 30 ml (for Butyrometers)	1 400 rpm	350 Xg
2	No. 2 Swing-out Rotor	4 x 50 ml (for BABCOCK TEST BOTTLES)	3 300 rpm	1 920 Xg



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com