



LABORATORY HIGH SPEED CENTRIFUGE BFD1BG1 (BCBH-311)

LABORATORY HIGH SPEED CENTRIFUGE BFD1BG1

HEMATOCRIT CENTRIFUGE



SPECIFICATIONS

Model	BFD1BG1
Old Model	BCBH-311
Max. RCF	15 800 Xg
Max. Speed	12 000 rpm
Rotor Capacity	24 Capillary tubes
Running Time	1 min - 99 min : 59 s
Drive System	Direct, brushless motor drive
Speed Accuracy	± 20 r/min
Noise Level	≤ 55 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	500 W
Dimension (W x D x H)	10.6 x 14.2 x 8.3 in / 27 x 36 x 21 cm
Weight with Rotor	40 lb / 14 kg
Alt Name	Hematocrit Centrifuge



ITEMS INCLUDED

No	Name
1	Rotor



1



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com