



LABORATORY LOW SPEED CENTRIFUGE BFD1BI1 (BCBLS-106)

LABORATORY LOW SPEED CENTRIFUGE BFD1BI1

BLOOD GROUP SEROLOGY CENTRIFUGE




SPECIFICATIONS

Model	BFD1BI1
Old Model	BCBLS-106
Max. RCF	1 790 Xg
Max. Speed	3 840 rpm
Max. Capacity	24 gel cards
Drive System	Direct, induction motor drive
Noise Level	≤ 50 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	100 W
Dimension (W x D x H)	14.8 x 14.2 x 11.8 in / 37.5 x 36 x 30 cm
Weight w/o Rotor	46 lb / 21 kg
Alt Name	Gel Card Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Rotor Capacity
1	No. 1 Swing-out Rotor	12 gel cards
2	No. 2 Swing-out Rotor	24 gel cards



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com